

Papain Dissociation of Cortical Neurons

McConnell Lab Protocol

Updated 29 May 2014

For use with the Worthington kit (<http://www.worthington-biochem.com/PDS/cat.html>)

Day Before:

- 1) Make ovomucoid-inhibitor solution ("High Inhibitory" or "HI"):
 - Add 32ml EBSS (vial 1) to albumin-ovomucoid mixture (vial 4)
 - Mix and equilibrate in incubator
 - Store at 4 deg.
 - Good for at least a few days
- 2) UV-treat hood overnight (alternatively, treat for one hour morning of FACS)
- 3) Make CSMN buffer (recipe below) (good for 2 days stored at 4 deg)
- 4) Treat plates with poly-D lysine
 - Coat 3 coverslips for each of 2x24-well plates with PDL (stock is 100X, add 100ul to 10ml water)
 - Let sit for 1 hour
 - Remove PDL and wash 2x with sterile H₂O
 - Add CSMN buffer to wells and store in incubator

Morning of:

- 1) Thaw Kynurenic acid (KYNA)
- 2) Place EBSS (provided in the kit) in the incubator
- 3) Remove embryos from pregnant dam and place in cold HBSS on ice.
Optional: Check for fluorescence under microscope.
- 4) Dissect embryos in cold HBSS. Store neocortical tissue in 15-ml Falcon tube on ice.
- 5) Place HI in the incubator with the lid slightly open for no longer than 35 minutes. (Note: do this just after dissection is complete.)

Dissociation:

- 1) Add 5-10ml EBSS to papain vial (vial 2). Place at 37 deg for 10 min. (Do this just before finishing the last dissection so papain has time to go into solution.)
- 2) Add 500ul EBSS to DNase vial (vial 3). Mix gently.
- 3) Make dissociation solution (Recipe for 1 dissociation):
 - Take 5ml papain in EBSS and add to new Falcon tube
 - Add 250ul DNase in EBSS
 - Add 5ul KYNA
- 4) Add all of the dissociation solution to tissue
- 5) Leave in incubator for 25-30 minutes, disturbing the tube every 5-10 minutes.
- 6) Gently triturate cells until they appear mostly dissociated, then centrifuge cells for 10 mins at #4 setting on McConnell centrifuge (~500g).

- 7) While centrifuging, make the low ovomucoid inhibitor solution ("Low Inhibitory" or "LI") (Recipe for 1 sample):
 - 2.7ml EBSS
 - 300ul HI
 - 150ul DNase
 - 3ul KYNA
- 8) Aspirate supernatant from cells.
- 9) Add all of the LI solution to the cells and carefully re-suspend (~3ml LI per sample)
- 10) In a new tube, add 5ml HI, and layer all of the cells in LI on top, for 8ml total
- 11) Centrifuge on #4 for ~7 mins
- 12) While centrifuging, add 20ul KYNA to 20ml CSMN buffer (or 40ul KYNA to 40ml CSMN).
- 13) Remove supernatant from cells.
- 14) Re-suspend cell pellet in ~500ul CSMN buffer and transfer to filter tube.
(Note: if >6 embryos were used, consider re-suspending in >500ul CSMN)
- 15) Spin filter cells using filter tube for 1 minute on #4. (Note: Each filter tube only holds 500ul, so may need to split sample into two tubes.)
- 16) Re-suspend filtered cells (already in 500ul+ CSMN buffer).
- 17) Spot 1ul onto each pre-FACS coverslip, wait 3 mins for cells to settle, then add 1-1.5 ml CSMN buffer, incubate at 37 deg
- 18) Sort in FACS facility.
Note: Cells seem to like being collected in a bacterial dish taped to the lid of 6-well plate containing 2-3 ml CSMN buffer (If more solution needs to be added to the cells being sorted, add CSMN **not** saline). If doing RNA or DNA analysis, collect cells in a microcentrifuge tube containing 500ul CSMN buffer or RNALater. Store cells on ice.
- 19) If culturing cells, collect everything from the dish into a falcon tube (~8ml total).

To Process After FACS:

- 1) Spin cells for 10 mins at #4 setting.
- 2) Carefully remove supernatant.
- 3) For cells in RNALater, resuspend in fresh RNALater and store at 4 deg.
- 4) For cells that were in CSMN, spot 1ul onto post-FACS coverslips, let sit 3 min, add 1-1.5ml CSMN buffer.
- 5) Incubate coverslips at 37 deg for 90 mins.
- 6) Remove CSMN buffer from plated cells and fix in 4% PFA for 5 mins.
- 7) Remove PFA and store in PBS at 4 deg.

CSMN buffer recipe (20 ml):

Ingredient	Amount	Final Concentration
DMEM	9.5 ml	50%
Neurobasal	9.5 ml	50%
Pyruvate	200 ul	1mM
L-Glutamine	200 ul	2mM
Insulin	200 ul	5ug/ml
Pen/Strep	200 ul	100U/100ug/ml
Sato	200 ul	1X
2M Glucose	350 ul	35mM
4% BSA	17 ul	0.34%